

microzone

microLYSIS[®] - Plus

For tough cells – e.g plants and tissues

Go from cells to PCR in less than 30 mins!

Protocol

Mix cells* with 20 µl microLYSIS-Plus

Overlay with mineral oil if necessary

Place in a Thermal Cycler

Lysis profile 1 - For very tough cells

Step 1: 65°C for 15 mins

Step 2: 96°C for 2 mins

Step 3: 65°C for 4 mins

Step 4: 96°C for 1 mins

Step 5: 65°C for 1 mins

Step 6: 96°C for 30 secs

Step 7: 20°C hold

A genuine Microzone product, made in the UK

Lysis profile 2** - For less tough cells

Step 1: 75°C for 5 mins

Step 2: 95°C for 2 mins

Step 3: 20°C hold

After lysis, all of the microLYSIS-Plus/DNA mixture can be used directly in PCR. It can make up to 40% of most PCR mixtures. Alternatively, it can be stored at -20°C for future use

*Free from PCR inhibitors

**Lysis profile 2 is not recommended if the released DNA is to be stored long term

Note: microLYSIS-Plus is a complex solution which releases (does not purify) the DNA from the cells. Estimation of DNA yield after lysis can therefore not be done using a spectrophotometer

Related Products

MegaMixes = Clear or blue PCR mixes for instant and accurate PCRs

microCLEAN = The easy way to clean-up your PCR fragments

BetterBase = Enhance/dilution buffer for sequencing

Supplied as a clear solution in multiple of 1 ml aliquots

Store at +4°C short term, -20°C long term

For Research Only